

CODAN

Infusion Pumps Docking Stations



Spare Parts
Accessories
Consumables



We care for every drop.



en

Preface

History

Document name	Version	Date	Main changes
15.597_A_ACCESS_en	A	2013-10-07	Initial version
15.597_B_CAT_en_Accessories	B	January 2016	Actualization, Harmonization, CODAN Corporate Identity

Edition notice

CODAN Infusionion pumps, Docking stations – Spare parts, Accessories, Consumables

The catalogue is for customers of the CODAN ARGUS AG and all those interested, who want to be informed about spare parts, accessories or consumables.

Intended use

The catalogue contains all orderable parts, accessories and consumables. They are listed with their complete name and reference number (important for orders). Most of them are illustrated by a picture.

If you want to place an order, please contact our sales people by e-mail or phone.

Contact address



CODAN ARGUS AG
Oberneuhofstrasse 10, 6340 Baar
Switzerland

Phone +41 41 785 09 44

Fax +41 41 785 09 40

www.codanargus.com

Support Technical Service Department
support@codanargus.com

Support Sales
codan@codanargus.com

Copyright

© Copyright: CODAN ARGUS AG. All rights reserved.


Terms and abbreviations

Item	Definition
A300 M&P, A500 M&P, A600 M&P	CODAN A300, A500, A600 M&P – Docking Stations
A60 M&P, A100 M&P	CODAN A60, A100 M&P – Docking Stations
A600S	CODAN A600S – Syringe Pump
A606S	CODAN A606S – Syringe Pump
A707V, A708V	CODAN A707V, A708V – Volumetric Infusion Pump
A717V, A718V	CODAN A717V, A718V – Volumetric Infusion Pump
AFF clamp	Anti-Free Flow Clamp
SSC	Safety Standard Check

Table of contents

Preface	2
Table of contents	4
Accessories	5
Consumables	12
Calibration Tools	15
CODAN A600S – Syringe Pump	17
CODAN A606S – Syringe Pump	23
CODAN A707V, A708V – Volumetric Infusion Pump	29
CODAN A717V, A718V – Volumetric Infusion Pump	36
CODAN A60, A100 M&P – Docking Station	47
CODAN A300, A500, A600 M&P – Docking Station	51

Accessories





Infusion stand, premium

- For single infusion pumps or docking stations
- Twin pole infusion stand
- 5 antistatic fork castors (2 with brakes)
- 4 bottle holders for infusions up to max. 1,5 kg each
- Mounting a docking station of max. 35 kg
- Mounting another device of max. 4 kg on twin pole

REF.: 601351



Infusion stand, standard

- For single infusion pumps
- Single pole infusion stand
- 5 antistatic double castors (2 with brakes)
- 4 bottle holders for infusions up to max. 1,5 kg each
- Mounting pumps of max. 5 kg

REF.: 601545



Bottle holder

- For CODAN volumetric pumps
- Available in 2 heights: 450 and 650 mm

REF.: 600987 / 450 mm

REF.: 601120 / 600 mm



Carrying handle

- For CODAN syringe pumps
- Pump remains compatible with docking station

REF.: 601235



Barcode scanner

- Cordless barcode imaging for healthcare applications
- Bluetooth compatible
- 1D and 2D barcodes are supported
- Protection class IP43

REF.: 601554

**Power supply for barcode scanner**

- PWRS-14000-058R 100-250VAC 5V/2A
- Delivered without power cable (please order the power cable separately according to your specific country needs)

REF.: 601564

**Power distributor**

- With 4 cables extendable to 1.80 m
- Standard plugs for medical devices
- With combination clamp for poles up to Ø 38 mm

REF.: 600450

**CODAN docking stations**

- For CODAN infusion pumps
- Available in 3 sizes:
- With 3 slots – A300 M&P
 - With 5 slots – A500 M&P
 - With 6 slots – A600 M&P
- Protection class II, IP21
 - Max. 100 VA
 - SSC after 24 months or 10 000 operating hours

Version P

- Power supply, 100 - 240 V

Version M

- Power supply, 100 - 240 V
- RS-232 Interface, Ethernet
- Interface for nurse call system

REF.: 601293 / A300P

REF.: 601294 / A500P

REF.: 601295 / A600P

REF.: 601296 / A300M

REF.: 601297 / A500M

REF.: 601298 / A600M

**Combination clamp**

- For CODAN infusion pumps
- For pole and rail fixation

REF.: 600978 / standard

REF.: 600979 / with rail holder V25x10

REF.: 600981 / with rail holder V35x8

REF.: 600982 / with T-rail holder V35x8

REF.: 600980 / with T-rail holder V35x10



Combination clamp for docking stations

- For pole and rail fixation

REF.: 601553 / A300, A500, A600 M&P



Rail holder for docking stations, V25x10

REF 601012 / A60 M&P, A100 M&P



Upgrade Kit

- Mountable onto combination clamp standard

REF.: 600983 / rail holder V25x10

REF.: 601223 / rail holder V35x8

REF.: 601224 / T-rails V35x8

REF.: 601225 / T-rails V35x10



Drop detector

- Component of CODAN volumetric pumps

REF.: 600986



Edgeboard hood

REF.: 600491



Cable RJ45 interface

- For CODAN infusion pumps
- For PC connection
- Update, configuration, read-out
- Software packages on demand

REF.: 600100



Adapter D-Sub for cable RJ45 interface

REF.: 600099

**Cable USB interface**

- For CODAN infusion pumps
- Updates, configurations, read-outs
- Software packages on demand

REF.: 601370

**Cable USB interface**

- Software packages on demand

REF.: 601496 / A300M, A500M, A600M

**Cable RS-232 interface**

- Cable length 3 m

REF.: 600586 / A60M, A100M

**Cable RS-242 cascade**

- Cable length 2 m

REF.: 600587 / A60M, A100M

**Cable nurse call**

- Cable length 2 m

REF.: 600985 / A707V, A708V

**Cable nurse call**

- Cable length 2 m

REF.: 600953 / A600S

**Cable nurse call**

- Cable length 2 m

REF.: 601483 / A717V, A718V, A606S,
A300M, A500M, A600M**Cable for external DC supply**

- For CODAN infusion pumps
- 20 V
- Cable length 3 m

REF.: 600989



Power cable US

- Connector 90° angled
- Cable length 2.5 m

REF.: 600576



Power cable Schuko

- Connector 90° angled
- Cable length 2.5 m

REF.: 600577



Power cable CH

- Cable length 2.5 m

REF.: 600578 / connector 90° angled

REF.: 600579 / connector straight



Power cable GB

- Connector 90° angled
- Cable length 2.5 m

REF.: 600575

REF.	Accessory If not marked different respective part is suitable to all CODAN products listed here.	A600S	A606S	A707V	A708V	A717V	A718V	A60M	A60P	A100M	A100P	A300M	A300P	A500M	A500P	A600M	A600P
600099	Adapter D-Sub for cable RJ45 interface	X	X	X	X	X	X										
600100	Cable RJ45 interface	X	X	X	X	X	X										
600450	Power distributor																
600491	Edgeboard hood	X	X	X	X	X	X										
600575	Power cable GB, 90°, 2.5 m																
600576	Power cable US, 90°, 2.5 m																
600577	Power cable Schuko, 90°, 2.5 m																
600578	Power cable CH, 90°, 2.5 m																
600579	Power cable CH, 90°, 2.5 m																
600586	Cable RS-232 interface, 3 m							X		X							
600587	Cable RS-242 cascade, 2 m							X		X							
600953	Cable nurse call, 2 m	X															
600978	Combination clamp standard	X	X	X	X	X	X	X	X	X	X						
600979	Combination clamp with rail holder V25x10	X	X	X	X	X	X										
600980	Combination clamp with T-rail holder V35x10	X	X	X	X	X	X										
600981	Combination clamp with rail holder V35x8	X	X	X	X	X	X										
600982	Combination clamp with T-rail holder V35x8	X	X	X	X	X	X										
600983	Upgrade kit: rail holder V25x10	X	X	X	X	X	X										
600985	Cable nurse call, 2 m			X	X												
600986	Drop detector			X	X	X	X										
600987	Bottle holder, 450 mm			X	X	X	X										
600989	Cable for external DC supply	X	X	X	X	X	X										
601012	Rail holder V25x10 for docking stations							X	X	X	X	X	X	X	X	X	X
601120	Bottle holder, 600 mm			X	X	X	X										
601223	Upgrade kit: for rail holder V35x8	X	X	X	X	X	X										
601224	Upgrade kit: for T-rails V35x8	X	X	X	X	X	X										
601225	Upgrade kit: for T-rails V35x10	X	X	X	X	X	X										
601235	Carrying handle	X	X														
601293	CODAN A300P - Docking Station																
601294	CODAN A500P - Docking Station																
601295	CODAN A600P - Docking Station																
601296	CODAN A300M - Docking Station																
601297	CODAN A500M - Docking Station																
601298	CODAN A600M - Docking Station																
601351	Infusion stand premium																
601370	Cable USB interface	X	X	X	X	X	X										
601483	Cable nurse call, 2 m		X			X	X					X		X		X	
601496	Cable USB interface for docking station											X		X		X	
601545	Infusion stand standard																
601553	Combination clamp for docking stations											X	X	X	X	X	X
601554	Barcode scanner					X	X										
601564	Power supply for barcode scanner					X	X										

Consumables



**Infusion set L86-P**

REF.: 601173 / A707V, A717V

**Infusion set LIGHT-SAVE™ L86-P**

REF.: 601513 / A707V, A717V

**Infusion set IVIP® 708/718 - V86**

REF.: 601217 / A708V, A718V

**Infusion set B88**

- For blood transfusion

REF.: 601176 / A707V, A717V

**Syringe CODAN 50 / 60 ml**

REF.: 601177

**Syringe CODAN 20 ml**

REF.: 601460

**Syringe CODAN 10 ml**


REF.: 601461

**Extension set E87-P, NoDEHP**

REF.: 601265

REF.	Consumable If not marked different respective part is suitable to all CODAN products listed here.						
		A600S	A606S	A707V	A708V	A717V	A718V
601173	Infusion set L86-P			X		X	
601176	Infusion set B88 for blood transfusion			X		X	
601177	Syringe CODAN 50 / 60 ml	X	X				
601217	Infusion set IVIP® 708/718 - V86				X		X
601265	Extension set E87-P, NoDEHP			X	X	X	X
601460	Syringe CODAN 20 ml	X	X				
601461	Syringe CODAN 10 ml	X	X				
601513	LIGHT-SAVE™ L86-P			X		X	

Calibration Tools





Pressure gauge

- 2,5 bar, with stopcock

REF.: 601506



Syringe calibration parts

- Part 1: Ø 31 mm, length 20 mm
- Part 3B: Ø 17 mm, length 120 mm

REF.: 601641 / part 1, A606S

REF.: 601639 / part 3B, A606S



Spring gauge for syringe pumps

REF.: 600904



Recommended balance for calibration

- Required scaling: 0.1 g

NOT AVAILABLE THROUGH CODAN

REF.	Calibration tool If not marked different respective part is suitable to all CODAN products listed here.						
		A600S	A606S	A707V	A708V	A717V	A718V
600904	Spring gauge for syringe pumps	X	X				
601506	Pressure gauge, 2,5 bar, with stopcock			X	X	X	X
601639	Calibration part 3B		X				
601641	Calibration part 1		X				



Spare Parts

CODAN A600S – Syringe Pump







Figure 1 A600S - Opened view, complete



Casing base V3

REF.: 600920




Casing cover

REF.: 600912

**Plastic foil panel**

REF.: 600770

**Display board**

REF.: 600938

**Cover driving head**

REF.: 600936

**Drive unit V3**

REF.: 600932

**Kit: flex cable V3**

REF.: 601164

**Side wall with DMS**


REF.: 600758

**Syringe clamp**

REF.: 600750

**Kit: pushbutton switch**

REF.: 600919



Main board

REF.: 600942 / V3.01

REF.: 600943 / V4.32



Power board

REF.: 600939 / 230 V

REF.: 600940 / 115 V



Battery cover 1 300 mAh

REF.: 600659



Battery cover 4 500 mAh

REF.: 600724



Lead-acid battery 6 V / 1 300 mAh

REF.: 600669



Lead-acid battery 6 V / 4 500 mAh

REF.: 600670



Syringe guide

REF.: 600744

**Female receptacle 3-pole**

REF.: 600915

**Case / male plug 7-pole**

REF.: 601186

**Motor complete**

REF.: 600922

**Seal for combination clamp**

REF.: 600705

**Fixing plate**

REF.: 600710

**Driving head V3**

- Without cover, delivered with an unmounted lid

REF.: 601620

**Holder with edgeboard**

REF.: 601015

**Working lever**

REF.: 600738

SPARE PARTS

CODAN A600S – Syringe Pump

REF.	Spare part
600659	Battery cover, 1 300 mAh
600669	Lead-acid battery 6 V / 1 300 mAh
600670	Lead-acid battery 6 V / 4 500 mAh
600705	Seal for combination clamp
600710	Fixing plate
600724	Battery cover, 4 500 mAh
600738	Working lever
600744	Syringe guide
600750	Syringe clamp
600754	Casing cover
600758	Side wall with DMS
600770	Plastic foil panel
600771	Label short instructions Danish
600772	Label short instructions German
600773	Label short instructions English
600774	Label short instructions French
600775	Label short instructions Dutch
600776	Label short instructions Italian
600777	Label short instructions Portuguese
600778	Label short instructions Swedish
600779	Label short instructions Spanish
600780	Label short instructions Czech
600781	Label plug external
600784	Label identification plate
600912	Casing cover
600915	Female receptacle 3-pole
600919	Kit: pushbutton switch
600920	Casing base V3
600922	Motor complete
600932	Drive unit V3
600936	Cover driving head
600938	Display board
600939	Power board, 230 V
600940	Power board, 115 V
600942	Main board V3.01
600943	Main board V4.32
601015	Holder with edgeboard
601164	Kit: flex cable V3
601186	Case / male plug 7-pole
601620	Driving head V3
601747	Kit: torsion spring (5 pieces)

Spare Parts

CODAN A606S – Syringe Pump






Figure 2 A606S - opened view, complete



Cover driving head

REF.: 600936



Driving head V3

- Without cover, delivered with an unmounted lid

REF.: 601620

**Kit: drive complete**

- Inclusive kit: fixation plate


REF.: 601707

**Main board**

REF.: 601271

**Casing cover**

REF.: 601060


**Casing base**

REF.: 601061 / 230 V

REF.: 601147 / 115 V

**Working lever**

REF.: 600738

**Kit: side wall with DMS**

- Inclusive kit: fixation plate

REF.: 601706

**Holder with edgeboard**

REF.: 601064

**Motor**

REF.: 601066

**Battery**

- 12 V / 1 800 mAh, NiMH

REF.: 601074

**Interface board RS-232**

REF.: 600555

**Kit: socket head cap screw M3x50 (10 pieces)**

REF.: 601740

**Kit: socket head cap screw M4x20 (10 pieces)**

REF.: 601741

**Kit: socket head cap screw M5x20 for combination clamp (10 pieces)**

REF.: 601744

**Kit: flex cable V3**

REF.: 601164

**Syringe clamp**

REF.: 600750

**Mains receptacle**

REF.: 600486

**Seal for combination clamp**

REF.: 600705

**Syringe guide**

REF.: 600744


**Plastic foil keyboard**

REF.: 600722

REF.	Spare part
600486	Mains receptacle
600555	Interface board RS-232
600658	Lithium coin CR1632.IB 3 V / 125 mAh (for main board)
600701	Battery cover
600705	Seal for combination clamp
600722	Plastic foil keyboard
600738	Working lever
600744	Syringe guide
600750	Syringe clamp
600787	Label short instructions Czech
600788	Label short instructions Swedish
600789	Label short instructions Portuguese
600790	Label short instructions Dutch
600791	Label short instructions Italian
600792	Label short instructions Danish
600793	Label short instructions English
600794	Label short instructions German
600881	Label identification plate 230 V
600882	Label identification plate 115 V
600883	Label short instructions Polish
600936	Cover driving head
600937	Driving head
601060	Casing cover
601061	Casing base 230 V
601064	Holder with edgeboard
601066	Motor
601074	Battery, 12 V / 1 800 mAh, NiMH
601103	Transformer 115 V
601104	Label short instructions French
601147	Casing base 115 V
601151	Label short instructions Turkish
601164	Kit: flex cable V3
601172	Label short instructions Slovak
601271	Main board
601620	Driving head V3
601706	Kit: side wall with DMS, inclusive: screw, washer, fixing plate
601707	Kit: drive complete
601735	Kit: pan head screw M3x10 (10 pieces)
601740	Kit: socket head cap screw M3x50 (10 pieces)
601741	Kit: socket head cap screw M4x20 (10 pieces)
601744	Kit: socket head cap screw M5x20 for combination clamp (10 pieces)
601745	Kit: fuse FST 125mA 5x20, 230V (10 pieces)
601746	Kit: fuse FST 250mA 5x20, 115V (10 pieces)
601747	Kit: torsion spring (5 pieces)
601748	Kit: non-slip rubber foot Ø10x3 black (16 pieces)

Spare Parts


CODAN A707V, A708V – Volumetric Infusion Pump





Display board

REF.: 600968



Sensor board

REF.: 600970 / A707V

REF.: 601052 / A708V



Main board

REF.: 601226 / A707V

REF.: 601234 / A708V



Combination clamp fixing element

REF.: 600132



Kit: socket head cap screw M4x35 (10 pieces)

REF.: 601742



Kit: socket head cap screw M5x20 for combination clamp (10 pieces)

REF.: 601744



Buzzer OBO

REF.: 600367

**Membrane for peristaltic**

REF.: 601732

**Cord sealing 2,0**

REF.: 600606

**Alarm window, right**

REF.: 600607

**Alarm window, left**

REF.: 600608

**Transformer**

REF.: 600956 / 230 V

REF.: 600961 / 115 V

**Insertion module AFF clamp**

REF.: 601057 / A708V

**Air detector**

REF.: 600976

**Cable stepper motor**

REF.: 601529



Pump unit

REF.: 600964 / A707V

REF.: 601050 / A708V



AFF clamp

REF.: 600966 / A707V



AFF board

REF.: 601157 / A708V



Pressure sensor Z1.0

REF.: 601188



Battery

- 12 V / 1 800 mAh, NiMH

REF.: 601259



Plastic foil panel

REF.: 600797 / A707V

REF.: 600815 / A708V



Edgeboard


REF.: 600969



Kit: non-slip rubber feet

- 16 x small, 16 x large

REF.: 601229

**Casing backside**

REF.: 600957 / 230 V

REF.: 600958 / 115 V

**Casing backside, blank**

REF.: 601230

**Casing frontside**

REF.: 601228

**Door**

REF.: 600960 / A707V

REF.: 601049 / A708V

**Mains receptacle**

REF.: 600486

SPARE PARTS

CODAN A707V, A708V – Volumetric Infusion Pump

REF.	Spare part If not marked different respective part is suitable to all CODAN products listed here.	A707V	A708V	
600132	Combination clamp fixing element			
600367	Buzzer OBO			
600486	Mains receptacle			
600606	Cord sealing, length 2.0 m			
600607	Alarm window, right			
600608	Alarm window, left			
600795	Label flow direction, 2-piece			
600797	Plastic foil panel	X		
600802	Label short instructions French			
600803	Label short instructions Italian			
600804	Label short instructions Czech			
600805	Label short instructions Portuguese			
600806	Label short instructions Swedish			
600807	Label identification plate	X		
600808	Label short instructions Dutch			
600809	Label short instructions Norwegian			
600810	Label short instructions Spanish			
600811	Label parking position drop detector			
600812	Label short instructions Danish			
600813	Label short instructions German			
600814	Label short instructions English			
600815	Plastic foil panel		X	
600818	Label identification plate		X	
600956	Transformer 230 V			
600957	Casing backside 230 V			
600958	Casing backside 115 V			
600960	Door			
600961	Transformer 115 V			
600964	Pump unit	X		
600966	AFF clamp			
600968	Display board			
600969	Edgeboard			
600970	Sensor board	X		
600976	Air detector			
601047	Casing frontside	X		
601049	Door	X		
601050	Pump unit		X	
601052	Sensor board		X	
601057	Insertion module AFF clamp	X		
601152	Label short instructions Turkish			
601157	AFF board			
601188	Pressure sensor Z1.0		X	
601226	Main board		X	
601228	Casing frontside		X	
601229	Kit: non-slip rubber feet (16 x small, 16 x large)			

REF.	Spare part If not marked different respective part is suitable to all CODAN products listed here.	A707V	A708V	
601230	Casing backside, blank			
601234	Main board		X	
601259	Battery, 12 V / 1.800 mAh, NiMH			
601529	Cable stepper motor	X		
601732	Membrane for peristaltic			
601742	Kit: socket head cap screw M4x35 (10 pieces)			
601743	Kit: socket head cap screw M4x65 A70X (10 pieces)			
601744	Kit: socket head cap screw M5x20 for combination clamp (10 pieces)		X	
601745	Kit: fuse FST 125mA 5x20, 230V (10 pieces)			
601746	Kit: fuse FST 250mA 5x20, 115V (10 pieces)			

Spare Parts

CODAN A717V, A718V – Volumetric Infusion Pump





Figure 3 Casing front, outside




Figure 4 Casing front, inside




Figure 5 Casing backside, outside




Figure 6 Casing backside, inside




Figure 7 Exploded view, complete

SPARE PARTS

CODAN A717V, A718V – Volumetric Infusion Pump

No.	REF.	Spare part If not marked different respective part is suitable to all CODAN products listed here.		
			A717V	A718V
1	601281	Casing frontside, complete	X	
	601485	Casing frontside, complete		X
2	601279	Plastic foil keyboard	X	
	601340	Plastic foil keyboard		X
3	601238	Alarm window, left		
4	601239	Alarm window, right		
5	600960	Door	X	
	601049	Door		X
7	601438	Main board with LCD	X	
	601439	Main board with LCD		X
8	600658	Lithium coin CR1632.IB 3 V / 125 mAh (for main board)		
9	600976	Air detector		
10	601188	Pressure sensor Z1.0		
11	601144	Insertion module AFF clamp		X
	600966	AFF clamp	X	
12	601491	Buzzer		
13	601615	Sensor board		
14	601343	Flat ribbon cable 6-pole		
15	601529	Cable stepper motor		
16	601275	Pump unit		
17	601732	Membrane for peristaltic		
18	601389	Casing backside 115 V		
	601274	Casing backside 230 V		
19	601282	Edgeboard		
20	601277	Power board		
21	601103	Transformer 115 V		
	601062	Transformer 230 V		
22	600486	Mains receptacle		
23	601748	Kit: non-slip rubber foot Ø10x3 black (16 pieces)		
24	601342	Flat ribbon cable 24-pole		
25	601074	Battery, 12 V / 1 800 mAh, NiMH		
26	601750	Kit: socket head cap screw M4x12 (12 pieces)		
27	601751	Kit: socket head cap screw M4x70 (12 pieces)		
28	600491	Edgeboard hood		
29	601744	Kit: socket head cap screw M5x20 for combination clamp (10 pieces)		
30	601749	Kit: screw cap (16 pieces)		
31	601745	Kit: fuse FST 125mA 5x20, 230V (10 pieces)		
	601746	Kit: fuse FST 250mA 5x20, 115V (10 pieces)		
32	600978	Combination clamp standard		

**Main board with LCD**

REF.: 601438 / A717V

REF.: 601439 / A718V

**Power board**

REF.: 601277

**AFF board**

REF.: 601157 / A718V

**Kit: socket head cap screw M5x20 for combination clamp (10 pieces)**

REF.: 601744

**Sensor board**

REF.: 601615

**Combination clamp fixing element**

REF.: 600132

**Kit: fuse (10 pieces)**

REF.: 601745 / Fuse FST 125 mA 5x20, 230 V

REF.: 601746 / Fuse FST 250 mA 5x20, 115 V

**Mains receptacle**

REF.: 600486

**Membrane for peristaltic**

REF.: 601732

**Kit: screw cap (16 pieces)**

REF.: 601749

**Air detector**

REF.: 600976

**Door**

REF.: 600960 / A717V

REF.: 601049 / A718V

**AFF clamp**

REF.: 600966 / A717V

**Insertion module AFF clamp**

REF.: 601144 / A718V

**Pressure sensor Z1.0**

REF.: 601188

**Alarm window, left**

REF.: 601238

**Alarm window, right**

REF.: 601239

**Flat ribbon cable 24-pole**

REF.: 601342

**Flat ribbon cable 6-pole**

REF.: 601343

**Pump unit**

REF.: 601275

**Cable stepper motor**

REF.: 601529

**Casing backside**

REF.: 601274 / 230 V

REF.: 601389 / 115 V



Casing frontside, complete

REF.: 601281 / A717V

REF.: 601485 / A718V



Plastic foil keyboard

REF.: 601279 / A717V

REF.: 601340 / A718V



Kit: socket head cap screw M4x12 (12 pieces)

REF.: 601750



Kit: socket head cap screw M4x70 (12 pieces)

REF.: 601751



Buzzer

REF.: 601491

REF.	Spare part If not marked different respective part is suitable to all CODAN products listed here.		
		A717V	A718V
600132	Combination clamp fixing element		
600486	Mains receptacle		
600658	Lithium coin CR1632.IB 3 V / 125 mAh (for main board)		
600795	Label flow direction, 2-piece		
600811	Label parking position drop detector		
600960	Door	X	
600966	AFF clamp		
600976	Air detector		
601049	Door		X
601062	Transformer 230 V		
601074	Battery, 12 V / 1 800 mAh, NiMH	X	
601103	Transformer 115 V	X	
601144	Insertion module AFF clamp		X
601157	AFF board		X
601188	Pressure sensor Z1.0		X
601238	Alarm window, left		
601239	Alarm window, right		
601274	Casing backside 230 V		
601275	Pump unit		X
601277	Power board		
601279	Plastic foil keyboard	X	
601281	Casing frontside, complete	X	
601282	Edgeboard		
601340	Plastic foil keyboard		X
601342	Flat ribbon cable 24-pole		
601343	Flat ribbon cable 6-pole		
601345	Label short instructions English	X	
601346	Label short instructions German	X	
601356	Label short instructions French		
601357	Label short instructions Italian		
601358	Label short instructions Danish		
601359	Label short instructions Dutch		X
601360	Label short instructions Norwegian		
601361	Label short instructions Swedish		
601363	Label short instructions Portuguese		
601364	Label short instructions Czech		
601365	Label short instructions Polish		
601366	Label short instructions Slovak		
601389	Casing backside 115 V		
601438	Main board with LCD	X	
601439	Main board with LCD		X
601477	Label identification plate 230 V	X	
601485	Casing frontside, complete		X
601491	Buzzer		
601514	Label identification plate 115 V	X	



SPARE PARTS

CODAN A717V, A718V – Volumetric Infusion Pump

REF.	Spare part If not marked different respective part is suitable to all CODAN products listed here.		
		A717V	A718V
601515	Label identification plate 230 V		X
601516	Label identification plate 115 V		X
601529	Cable stepper motor	X	
601615	Sensor board		X
601732	Membrane for peristaltic		
601744	Kit: socket head cap screw M5x20 for combination clamp (10 pieces)		X
601745	Kit: fuse FST 125mA 5x20, 230V (10 pieces)		
601746	Kit: fuse FST 250mA 5x20, 115V (10 pieces)	X	
601748	Kit: non-slip rubber foot Ø10x3 black (16 pieces)		X
601749	Kit: screw cap (16 pieces)		X
601750	Kit: socket head cap screw M4x12 (12 pieces)		
601751	Kit: socket head cap screw M4x70 (12 pieces)		

Spare Parts

CODAN A60, A100 M&P – Docking Station



CODAN A60, A100 M&P – Docking Station

**Main board**

REF.: 601005 / A60P, A100P

**Main board**


REF.: 601007 / A60M, A100M

**Gender changer DB9 m/m**

REF.: 600492

**Centering fixation for docking station**

REF.: 600497

**Cover 103**

REF.: 600498 / A60 M&P

**Locking device for docking station**

- Screw: headless screw M3x3

REF.: 600499

**Coil spring 1 for docking station**

REF.: 600511

**Transformer holder**

REF.: 600515

**Voltage regulator**

REF.: 600993



Handle with axle

REF.: 600538



Casing backside A

REF.: 600540 / A100 M&P



Casing backside 3

REF.: 600541 / A60 M&P



Gender changer D-Sub 9-pole F/F

REF.: 600898




Gender changer D-Sub 9-pole M/F

REF.: 600899



Null modem adapter D-Sub 9-pole M/F

REF.: 600900



Transformer

REF.: 600998 / 230 V

REF.: 600999 / 115 V



Cover backside, complete

REF.: 601001 / A100 M&P



Cover backside, complete

REF.: 601002 / A60 M&P


SPARE PARTS

CODAN A60, A100 M&P – Docking Station

REF.	Spare part If not marked different respective part is suitable to all CODAN products listed here.	A60M	A60P	A100M	A100P
600492	Gender changer DB9 m/m	X		X	
600497	Centering fixation for docking station			X	X
600498	Cover 103	X	X		
600499	Locking device for docking station (screw: headless screw M3x3)				
600511	Coil spring 1 for docking station				
600515	Transformer holder				
600538	Handle with axle				
600540	Casing backside A			X	X
600541	Casing backside 3	X	X		
600821	Label ARGUS			X	
600822	Label ARGUS				X
600823	Label ARGUS	X			
600824	Label ARGUS		X		
600825	Label identification plate			X	
600826	Label identification plate				X
600827	Label identification plate, 115 V				X
600828	Label identification plate	X			
600829	Label identification plate		X		
600830	Label identification plate 115 V		X		
600898	Gender changer D-Sub 9-pole F/F	X		X	
600899	Gender changer D-Sub 9-pole M/F	X		X	
600900	Null modem adapter D-Sub 9-pole M/F	X		X	
600993	Voltage regulator				
600998	Transformer 230 V				
600999	Transformer 115 V				
601001	Cover backside, complete			X	X
601002	Cover backside, complete	X	X		
601004	Extended board			X	X
601005	Main board		X		X
601007	Main board	X		X	

Spare Parts

CODAN A300, A500, A600 M&P – Docking Station






Figure 8 Exploded view 1 - M docking station




Figure 9 Exploded view 2 - M docking station




Figure 10 Exploded view - P docking station




Figure 11 Exploded views - details

SPARE PARTS

CODAN A300, A500, A600 M&P – Docking Station

No.	REF.	Spare part If not marked different respective part is suitable to all CODAN ARGUS products listed here.	A300M	A300P	A500M	A500P	A600M	A600P	
1	601744	Kit: socket head cap screw M5x20 for combination clamp (10 pieces)							
2	601553	Combination clamp for docking station							
3	601737	Kit: screw PT-LS KA40x10 (10 pieces)							
4	601311	Cover backside					X	X	
	601312	Cover backside			X	X			
	601313	Cover backside	X	X					
5	601384	Insulation foil							
6	601613	Main board (screws: 5 x PT-LS-KA25x8)							
7	601291	Connection ribbon cable – 20-pole							
8	601303	Distribution board with 2 slots (screws: 6 x PT-LS-KA25x8)			X	X			
	601304	Distribution board with 3 slots (screws: 6 x PT-LS-KA25x8)					X	X	
9	600479	Equipotential plug POAG 6/11							
10	601339	Cover for power supply (screws: 2 x PT-LSKA25x6)							
11	601288	Power supply 115 – 230 V							
12	601503	Mains socket – M docking stations	X		X		X		
	601612	Mains socket – P docking stations		X		X		X	
13	601308	Enclosure housing					X		
	601309	Enclosure housing			X				
	601310	Enclosure housing	X						
14	601302	Power board with 3 slots (screws: 6 x PT-LS-KA25x8)							
15	601531	Insulation foil 2							
16	601290	Connection ribbon cable – 40-pole							
17	600538	Handle with axle							
18	601337	LED light guide							
19	601479	Serial connector internal cable	X		X		X		
20	600499	Locking device for docking station (screw: headless screw M3x3)							
21	600511	Coil spring 1 for docking station							
22	601534	Solder lug 6.3 x 0.8							
23	601542	Functional earth cable – serial connector							
24	601543	Functional earth cable – mains connector							
25	601544	Functional earth cable – main board							
26	601734	Kit: toothed lock washer M6 (10 pieces)							
27	601739	Kit: hexagon nut M6 (10 pieces)							
28		Screw PT-LSKA25x8							
29	601338	Mounting frame for power supply							
30	601736	Kit: screw PT-LS KA30x8 (10 pieces)							
31	601292	Connection cable power cable							
32	601348	Top cover	X	X	X	X			
33	601505	Jumper board							
34	601386	Power cable LED							
35	601305	Enclosure housing						X	
	601306	Enclosure housing				X			
	601307	Enclosure housing		X					

**Main board**

REF.: 601613 / A300M, A500M, A600M

**Power board 3 slots**

REF.: 601302

**Distribution board 2 slots**

REF.: 601303 / A500 M&P

**Distribution board 3 slots**

REF.: 601304 / A600 M&P

**Kit: toothed lock washer M6 (10 pieces)**

REF.: 601734

**Kit: screw PT-LS KA30x8 (10 pieces)**

REF.: 601736

**Kit: screw PT-LS KA40x10 (10 pieces)**

REF.: 601737

**Kit: hexagon nut M6 (10 pieces)**

REF.: 601739

**Kit: socket head cap screw M5x20 for combination clamp (10 pieces)**

REF.: 601744

**Equipotential plug POAG 6/11**

REF.: 600479

**Locking device for docking station**

- Screw: headless screw M3x3

REF.: 600499

**Connection cable power cable**

REF.: 601292

**Coil spring 1 for docking station**

REF.: 600511

**Handle with axle**

REF.: 600538

**Power supply**

- 100 V – 240 V

REF.: 601288

**Connection ribbon cable 40-pole**

REF.: 601290

**Connection ribbon cable 20-pole**

REF.: 601291

**Enclosure housing**

REF.: 601305 / A600P

REF.: 601308 / A600M

**Enclosure housing**

REF.: 601306 / A500P

REF.: 601309 / A500M

**Enclosure housing**

REF.: 601307 / 300P

REF.: 601310 / 300M

**Cover backside**

REF.: 601311 / A600 M&P

**Cover backside**

REF.: 601312 / A500 M&P

**Cover backside**

REF.: 601313 / A300 M&P

**LED light guide**

REF.: 601337

**Mounting frame for power supply**

REF.: 601338

**Cover for switching power supply**

- Screws: 2 x PT-LSKA25x6

REF.: 601339

**Insulation foil**

REF.: 601384

**Top cover**

REF.: 601348 / A300, A500 M&P

**Power cable LED**

REF.: 601386

**Serial connector internal cable**

REF.: 601479 / A300M, A500M, A600M


**Solder lug**

- 6.3 x 0.8 mm


REF.: 601534

**Mains socket**

REF.: 601503 / A300M, A500M, A600M

**Jumper board**

REF.: 601505

**Insulation foil 2**

REF.: 601531

**Functional earth cable**

- Serial connector

REF.: 601542

**Functional earth cable**

- Mains connector

REF.: 601543

**Functional earth cable**

- Main board

REF.: 601544

**Mains socket**

REF.: 601612 / A300P, A500P, A600P

REF.	Spare part If not marked different respective part is suitable to all CODAN products listed here.	A300M	A300P	A500M	A500P	A600M	A600P
600479	Equipotential plug POAG 6/11						
600499	Locking device for docking station (screw: headless screw M3x3)						
600511	Coil spring 1 for docking station						
600538	Handle with axle						
601288	Power supply, 100 – 240 V						
601290	Connection ribbon cable – 40-pole						
601291	Connection ribbon cable – 20-pole						
601292	Connection cable power cable						
601302	Power board 3 slots (screws: 6 x PT-LS-KA25x8)						
601303	Distribution board 2 slots (screws: 6 x PT-LS-KA25x8)						
601304	Distribution board 3 slots (screws: 6 x PT-LS-KA25x8)						
601305	Enclosure housing					X	
601306	Enclosure housing			X			
601307	Enclosure housing		X				
601308	Enclosure housing				X		
601309	Enclosure housing			X			
601310	Enclosure housing	X					
601311	Cover backside				X	X	
601312	Cover backside			X	X		
601313	Cover backside	X	X				
601337	LED light guide			X			
601338	Mounting frame for power supply	X					
601339	Cover for power supply (screws: 2 x PT-LSKA25x6)				X	X	
601348	Top cover			X	X		
601384	Insulation foil	X	X				
601386	Power cable LED						
601479	Serial connector internal cable						
601503	Mains socket – M docking stations						
601505	Jumper board	X	X	X	X		
601531	Insulation foil 2						
601534	Solder lug 6.3 x 0.8						
601542	Functional earth cable – serial connector	X		X		X	
601543	Functional earth cable – mains connector	X		X		X	
601544	Functional earth cable – main board						
601612	Mains socket – P docking stations		X		X	X	
601613	Main board (screws: 5 x PT-LS-KA25x8)						
601734	Kit: toothed lock washer M6 (10 pieces)						
601736	Kit: screw PT-LS KA30x8 (10 pieces)						
601737	Kit: screw PT-LS KA40x10 (10 pieces)						
601739	Kit: hexagon nut M6 (10 pieces)		X		X	X	
601744	Kit: socket head cap screw M5x20 for combination clamp (10 pieces)						

CODAN Worldwide



CODAN is known internationally as a manufacturer and supplier of disposable medical transfer systems. The CODAN Companies have more than 1500 employees around the world.

The name CODAN is synonymous with reliability, quality and precision based on the know-how and experience gained from more than 50 years of research and development.

Company-owned production facilities and sales companies around the world are a guarantee for efficient production, a tight-knit sales network and a first-class service for customers in the healthcare sector.

CODAN Product range

- Infusion sets
- Transfusion sets
- Extension lines and manifold connectors
- Infusion and transfusion accessories
- Infusion filters and filter systems
- Neonatology/Paediatric products
- Withdrawal, preparation and administration systems
- CODAN Cyto®
- Chemoprotect® products
- CODAN ENTERAL
- Disposable syringes
- Invasive blood pressure monitoring systems
- Infusion pumps
- Continence care products
- Other CODAN Products

CODAN Companies

CODAN Medizinische Geräte GmbH & Co KG · Deutschland

CODAN pvb Critical Care GmbH · Deutschland

CODAN pvb Medical GmbH · Deutschland

CODAN PORTUGAL, S.A. · Portugal

CODAN 11, S.A. · Portugal

CODAN US Corporation · California · USA

CODAN Inc. · California · USA

CODAN NORGE AS · Norge

CODAN TRIPLUS AB · Sverige

CODAN Limited · Great Britain

CODAN FRANCE Sarl · France

CODAN Medical AG · Schweiz

CODAN ARGUS AG · Schweiz

CODAN BV · Nederland

CODAN s.r.l. · Italia

CODAN Medical GmbH · Österreich

CODAN Steritex ApS · Danmark

CODAN Medical ApS · Danmark

CODAN DEHA ApS · Danmark

Compliance of the established quality management systems with the provisions of EN ISO 13485 and the Council Directive 93/42/EEC has been certified by the relevant, competent notified bodies. Our products are thus entitled to bear the following marks:

CE 0123 TÜV Product Service
CODAN Medizinische Geräte GmbH & Co KG
D-23738 Lensahn

CE 0543 Presafe Denmark A/S
CODAN Steritex ApS · DK-3060 Espergærde

CE 0543 Presafe Denmark A/S
CODAN Medical ApS · DK-4970 Rødby

CE 0123 TÜV Product Service
CODAN pvb Critical Care GmbH
D-85661 Forstinning



CODAN ARGUS AG

Oberneuhofstrasse 10, 6340 Baar, Switzerland

Phone +41 41 785 09 44, Fax +41 41 785 09 40

codan@codanargus.com, www.codanargus.com

We care for every drop.

